

# Production Engineering - Medical Equipment Division

[PEMED HOME](http://www.pemed.com) 6035 E. 38th Avenue, Denver, Colorado 80207 USA (303) 393-7800 FAX (303) 393-1482

## Carts

### Endoscopy cart with monitor - \$675



List: \$4,200. MedCon: VG. Stryker Smart cart, model 240-86. 115 V, 60 Hz.

Features:

- Sony Trinitron 19 inch color video monitor, model PVM-1943MD

on rotating monitor stand with wrap-around restraint

- 2 locking Plexiglas front windows
- 3 stationary and one pull-out shelves each with back cord wrap
- internal vertical cord guides
- 2 AC power bars
- chart holder on side
- side door for equipment access
- bottom storage drawer
- side transport handles
- cord hangar in back
- front locking casters
- 13 foot heavy duty power cord with hospital grade plug



Overall dimensions of the unit:

- width: 34.5 in.
- depth: 31.5
- height: 53 plus 18 for monitor

Shipping weight is about 300 pounds.

---

## Supply cart - \$245



List: \$650. MedCon: VG. Knape-Vogt model 1330 cart.

Features:

- rolling stainless steel cart on four casters
- three 3.35 in. and one 7.75 in. deep drawers
- retaining ring on top work area

Overall dimensions of the unit:

- width: 18.75 in.
- depth: 15.75
- height: 32.5

Shipping weight is about 60 pounds.

---

## Endoscopy tower cart - \$675



List: \$3,100. MedCon: VG. Wolf model 31121.101 cart

Features:

- 4 rack shelves plus one pullout shelf
- padded equipment braces
- gas tank holder
- power cords channel

- back transport handle
- bottom drawer (locking)
- cord wrap

Overall dimensions of the unit:

- width: 25.5 in.
- depth: 28
- height: 58.75

Shipping weight is about 100 pounds.

---

## Supplies and records cart - \$575



List: \$980. MedCon: VG. Datel cart. Mobile workstation for patient records and clinic supplies.

Features:

- rolling cart on four casters with transport handles

on each side and bottom bumpers

- 36 x 18.75 in. top work surface with lip
- 2 full length locking top drawers, 2.75 in. deep
- 2 half length bottom drawers, 9.75 in. deep
- adjustable dividers included in drawers and records storage area
- power bar, six AC outlets, 125 V, 60 Hz.

Overall dimensions of the unit:

- width: 41.5 in.
- depth: 22
- height: 42.5

Shipping weight is about 180 pounds.

---

## Stainless steel cart - \$475



List: \$1,800. MedCon: VG. Pedigo Steril-Gard cart.

Features:

- top surface is 41 x 36 in.
- internal shelf is 21 x 36 in and can be mounted at eight positions

on the bracket strip

- open doors swing all the way to side of cart and latch in place
- two opposite corner casters have locking levers

Outside dimensions of the cart:

- width: 46 in.
- depth: 28.5
- height: 39

Interior dimensions of the cart:

- width: 40 in.
- depth: 24
- height: 29.5

Shipping weight is about 180 pounds.

## Locking pharmaceutical cart - \$475



List: \$2,900. MedCon: VG. locking pharmaceutical cart, with four compartments, two on each side. Compartments lock, although purchaser will need to re key the lock.

Cart is about 5 feet tall, 3.5 feet wide, 2 1/2 feet deep.

Shipping weight is about 300 pounds.

## Surgery storage cart - \$750



List: \$3,200. MedCon: ex. Surgical storage cart, with adjustable height wire rack shelves, and 18 supply storage shelves.

Appears to be made by either Hillrom or Metro.

Overall dimensions:

- height: 74"
- length: 57"
- working depth: 19"
- overall depth: 27"

Shipping weight is about 300 pounds, although the cart itself weighs in at about 175 pounds.

## Endoscopy cart, monitor and printer - \$1,275



List: \$6,000. MedCon: ex. Endoscopy cart.

Features:

- storage drawer
- 4 shelves for equipment and instruments
- monitor support as shown
- includes high resolution video monitor, video printer, and power bar

Shipping weight is about 300 pounds.

## Endoscopy cart - \$675



List: \$2,700. MedCon: VG. Wolf endoscopy cart, model 2191.201.

Features:

- keyed lift-in glass doors with magnetic clasps form two

compartments with two shelves (20,5 x 21.5 in.) each

- each compartment has two sets of three AC outlets forming a common power bar, 120 V, 50/60 Hz.
- 24.5 x 22.5 in. lipped top with port to accommodate monitor cords
- storage drawer at bottom
- keyed back access panel
- transport handles
- two locking casters
- cord wrap

Overall dimensions of the unit:

- width: 29.25 in.
- depth: 27
- height: 52.5

Shipping weight is about 200 pounds.

## Endoscopy surgery carts - \$675 each

### 1. Olympus



### 2. Stryker smart cart



3. Stryker, 1/2 height Plexiglas door, four shelves, monitor holder



4. Stryker, full height Plexiglas door, four shelves, monitor holder



Shipping weight is about 250 pounds each.

## Stainless steel shelf and cabinet - \$1,375 each



List: \$3,925. MedCon: EX. Stainless steel shelf with storage cabinet. All stainless steel construction. Made by Shampagne medical products. Unit has one work surface, two open shelves with rails, and two shelves in an enclosed locking cabinet. On castors.

Dimensions:

- width: 37 inches
- depth: 18.5 inches
- height: 64 inches

Shipping weight is about 175 pounds each. Three units in stock.


## Stainless steel surgical storage cabinet - \$1,675



List: \$7,300. MedCon: ex. stainless steel storage cabinet, with adjustable interior shelves, a sloped top and sliding glass doors. may be used for storage of disposables, blankets, IV solutions, wrapped sterile instruments or other ICU, OR and surgical items.

Overall outside dimensions:

- length: 24 inches
- width: 60 inches
- height: 86 inches



Inside dimensions, bottom compartment:

- length: 23"
- width: 59 3/4
- height: 86

Inside dimensions, top compartment:

- length: 23 inches
- width: 59 3/4
- height: 42 1/2

Shipping weight is about 400 lbs.

---

## Radiolucent lab cart - \$375



List: \$900. MedCon: VG. LabConco lab cart. With radiolucent top.

Masonite top, steel frame, castors.

Shipping weight is about 85 pounds.

---

## Stainless steel case carts - \$475 each

List: \$3,450 each. MedCon: VG. Stainless steel case carts.



Features:

- 2 interior full length adjustable height stainless shelves
- dual locking doors
- on 4 steerable castors
- Overall dimensions are about 4 feet wide, 4 feet high, 2 feet deep.

Shipping weight is about 300 pounds each. 6 units in stock.

---

## ESU cart - \$375 each



List: \$622. MedCon: G. Valleylab ESU cart.

3 units in stock.

Shipping weight is about 75 pounds.

---

## Case / supply cart - \$250



List: \$665. MedCon: EX. Manufacturer is unknown, but probably Metro. Molded plastic cart in excellent condition. With folding swing out door and several adjustable height shelves.

On 4 castors. 32" long, 22" deep, 42" high.

Shipping weight is about 100 pounds.

---

## Endoscopic video cart - \$450

List: \$1,900. MedCon: EX. 6' high overall.



Features:

- 8 power outlets for power distribution
  - doors on both front and back of unit
  - 4 locking castors
  - 2 support arms for equipment that hang off of top of unit
  - one of the support arms has adjustable height
  - 3 internal shelves
- 

© Copyright 1993 - 2009 by Production Engineering. PEMED® is a registered trademark of Production Engineering. MedCon is also a trademark of Production Engineering. All rights reserved. Contact PEMED for license requests. The usual copyright notices and disclaimers. Except for search engines or link listing, no part of this website, including derived information, may be republished without express written permission and payment of fees.

---

# Production Engineering - Medical Equipment Division

6035 E. 38th Avenue, Denver, Colorado 80207 USA (303) 393-7800 FAX (303) 393-1482

---

80207  
(303) 393-7800  
(303) 393-1482  
www.pemed.com  
markz@pemed.com

6035 E. 38th Avenue  
Denver, Colorado 80207  
USA  
telephone (303) 393-7800  
FAX (303) 393-1482  
<http://www.pemed.com>  
email: [markz@pemed.com](mailto:markz@pemed.com)

6035 E. 38th Avenue  
Denver, Colorado 80207  
USA  
telephone (303) 393-7800  
FAX (303) 393-1482  
<http://www.pemed.com>  
email: [markz@pemed.com](mailto:markz@pemed.com)