

Production Engineering - Medical Equipment Division

[PEMED HOME](#) 6035 E. 38th Avenue, Denver, Colorado 80207 USA (303) 393-7800 FAX (303) 393-1482

Carts, Cabinets and Stainless Steel

To get this entire section as a printable PDF document [click here](#)

Stainless steel cart - \$450



List: \$950. MedCon: VG. Blickman utility cart.

Features:


- stainless steel construction
- top shelf (35 x 21.25 inches) at 35.5 inches high with containing rim; bottom shelf at 14 inches high
- rubber bumper at bottom (shows wear)
- four castors

Overall dimensions of the unit:

- width: 42.5 in.
- depth: 24.75
- height: 37

Shipping weight is about 55 pounds.

Equipment or supplies cart - \$475



List: \$1,180. MedCon: VG. Stainless Equipment Company utility cart.

Features:

- stainless steel construction
- top shelf (23.5 x 38 inches) at 42 inches high
- middle and bottom shelves (20.5 x 36 inches) at 8 and 24 inches high
- transport handle
- four castors

Overall dimensions of the unit:

- width: 40 in.
- depth: 24.25
- height: 43.5

Shipping weight is about 70 pounds.

Stainless cabinet - \$1,750



List: \$6,272. MedCon: VG. Medical supplies cabinet, stainless steel construction.

Features:

- stainless steel construction
- upper cabinet with windowed door and two glass shelves, 28.5 x 14.5 x 13/32 inches
- padlock can be attached to upper cabinet
- bottom compartment with one stainless steel shelf
- quarter-turn handles

Overall dimensions of the unit:

- width: 30.5 in.
- depth: 17.75
- height: 59.5

Shipping weight is about 220 pounds.

Endoscopy cart - \$375



List: \$900. MedCon: VG. Olympus Keymed model TCH-1 endoscopy cart.

Features:


- four 16-inch wide shelves with lips or retaining bar
- extendable pole for hanging fluid bags
- supply bin
- two locking castors
- rear transport handle
- cord wrap
- bottom bumpers

Overall dimensions of the unit:

- width: 22.5 in.
- depth: 23
- height: 54

Shipping weight is about 80 pounds.

Procedure cart - \$675 each



List: \$1,819. MedCon: VG. HOMED treatment cabinet.

Features:

- locking cabinet with two drawers
- transport handle on left side
- optional retaining rail

Overall dimensions of the unit:

- width: 29.5 in.
- depth: 18.5
- height: 43.5 (plus 3 inches for retaining rail)

Shipping weight is about 112 pounds. 4 units in stock.

Motorized equipment cart - \$1,275



List: \$9,800. MedCon: VG. Basically this is a cart with 1 motorized axle. Varian Medical Systems Motorized equipment cart, ref. AL20030000. The cart appears to be designed for computers and printing for radiotherapy systems but should be adaptable for many uses.

Features:

- two 12V batteries in series provide 24VDC power source and charge from 100-240 V, 50-60 Hz, current; does not include charger.
- right top and left top openings (slide-back fold-up doors) for accessing equipment are 15.5 and 20 inches
- right front pull up panel provides 16 inch access
- bottom left area with hinged door could be used for a computer
- tilt bins in front and back can be used for patient records
- control panel has stop/reset button to prevent accidental movement, lighted forward/reverse direction indicator, and battery charge indicator
- pressure on handle activates movement speed



Overall dimensions of the unit:

- width: 55 in.
- depth: 31
- height: 43

Shipping weight is about 650 pounds.

Endoscopy cart with monitor - \$675



List: \$4,200. MedCon: VG. Stryker Smart cart, model 240-86. 115 V, 60 Hz.

Features:

- Sony Trinitron 19 inch color video monitor, model PVM-1943MD

on rotating monitor stand with wrap-around restraint

- 2 locking Plexiglas front windows
- 3 stationary and one pull-out shelves each with back cord wrap
- internal vertical cord guides
- 2 AC power bars
- chart holder on side
- side door for equipment access
- bottom storage drawer
- side transport handles
- cord hangar in back
- front locking casters
- 13 foot heavy duty power cord with hospital grade plug



Overall dimensions of the unit:

- width: 34.5 in.
- depth: 31.5
- height: 53 plus 18 for monitor

Shipping weight is about 300 pounds.

Stainless steel cart- \$475

List: \$1,800. MedCon: VG. Pedigo Steril-Gard cart.



Features:

- top surface is 41 x 36 in.
- internal shelf is 21 x 36 in and can be mounted at eight positions on the bracket strip
- open doors swing all the way to side of cart and latch in place
- two opposite corner casters have locking levers

Outside dimensions of the cart:

- width: 46 in.
- depth: 28.5
- height: 39

Interior dimensions of the cart:

- width: 40 in.
- depth: 24
- height: 29.5

Shipping weight is about 180 pounds.

Locking pharmaceutical cart - \$475



List: \$2,900. MedCon: VG. locking pharmaceutical cart, with four compartments, two on each side. Compartments lock, although purchaser will need to re key the lock.

Cart is about 5 feet tall, 3.5 feet wide, 2 1/2 feet deep.

Shipping weight is about 300 pounds.

Surgery storage cart - \$750



List: \$3,200. MedCon: ex. Surgical storage cart, with adjustable height wire rack shelves, and 18 supply storage shelves.

Appears to be made by either Hillrom or Metro.

Overall dimensions:

- height: 74"
- length: 57"
- working depth: 19"
- overall depth: 27"

Shipping weight is about 300 pounds, although the cart itself weighs in at about 175 pounds.

Endoscopy cart, monitor and printer - \$1,275



List: \$6,000. MedCon: ex. Endoscopy cart.

Features:

- storage drawer
- 4 shelves for equipment and instruments
- monitor support as shown
- includes high resolution video monitor, video printer, and power bar

Shipping weight is about 300 pounds.

Endoscopy cart - \$675



List: \$2,700. MedCon: VG. Wolf endoscopy cart, model 2191. 201.

Features:

- keyed lift-in glass doors with magnetic clasps form two

compartments with two shelves (20,5 x 21.5 in.) each

- each compartment has two sets of three AC outlets forming a common power bar, 120 V, 50/60 Hz.
- 24.5 x 22.5 in. lipped top with port to accommodate monitor cords
- storage drawer at bottom
- keyed back access panel
- transport handles
- two locking casters
- cord wrap

Overall dimensions of the unit:

- width: 29.25 in.
- depth: 27
- height: 52.5

Shipping weight is about 200 pounds.

Endoscopy surgery carts - \$675 each

1. Olympus



2. Stryker smart cart



3. Stryker, 1/2 height Plexiglas door, four shelves, monitor holder



4. Stryker, full height Plexiglas door, four shelves, monitor holder



Shipping weight is about 250 pounds each.

Radiolucent lab cart - \$375



List: \$900. MedCon: VG. LabConco lab cart. With radiolucent top.

Masonite top, steel frame, castors.

Shipping weight is about 85 pounds.

Stainless steel case carts - \$475 each

List: \$3,450 each. MedCon: VG. Stainless steel case carts.



Features:

- 2 interior full length adjustable height stainless shelves
- dual locking doors
- on 4 steerable castors
- Overall dimensions are about 4 feet wide, 4 feet high, 2 feet deep.

Shipping weight is about 300 pounds each. 6 units in stock.

ESU cart - \$375 each



List: \$622. MedCon: G. Valleylab ESU cart.

3 units in stock.

Shipping weight is about 75 pounds.

Case / supply cart - \$250



List: \$665. MedCon: EX. Manufacturer is unknown, but probably Metro. Molded plastic cart in excellent condition. With folding swing out door and several adjustable height shelves.

On 4 castors, 32" long, 22" deep, 42" high.

Shipping weight is about 100 pounds.

Endoscopic video cart - \$450

List: \$1,900. MedCon: EX. 6' high overall.



Features:

- 8 power outlets for power distribution
- doors on both front and back of unit
- 4 locking castors
- 2 support arms for equipment that hang off of top of unit
- one of the support arms has adjustable height
- 3 internal shelves

© Copyright 1993 - 2011 by Production Engineering. PEMED® is a registered trademark of Production Engineering. MedCon is also a trademark of Production Engineering. All rights reserved. Contact PEMED for license requests. The usual copyright notices and disclaimers. Except for search engines or link listing, no part of this website, including derived information, may be republished without express written permission and payment of fees.

Production Engineering - Medical Equipment Division

6035 E. 38th Avenue, Denver, Colorado 80207 USA (303) 393-7800 FAX (303) 393-1482