

Production Engineering - Medical Equipment Division

[PEMED HOME](#) 6035 E. 38th Avenue, Denver, Colorado 80207 USA (303) 393-7800 FAX (303) 393-1482

Physical Therapy / Rehab

To get this entire section as a printable PDF document [click here](#)

Washcloth warmer - \$395 each



List: \$1,900. MedCon: VG. Sage Products, Inc. model 7934 warmer for Sage Comfort Bath cleansing washcloths. 120 V, 60 Hz.

Features:

- 14 count capacity
- heats to 129°F (tested with empty chamber)
- orange LED digital temperature display
- lighted indicators for product ready or not ready
- power on/off Rocker switch with red indicator light
- 6-foot power cord with hospital grade plug

Overall dimensions of the unit:

- width: 14 in.
- depth: 14.75
- height: 17.5

Shipping weight is about 30 pounds. 3 units in stock.

Hydrotherapy bath - \$975



List: \$4,900. MedCon: VG. Arjo Sidekick dual tub hydrotherapy bath. 115 V, 60 Hz.

Features:

- unit can be wheeled to patients for hydro massage of upper or lower limbs
- bath height raises from 19 to 42 inches by electromechanical lift; pendant with raise and lower buttons
- each tub has water jets with aeration
- digital temperature display
- drain hose and power out pump for emptying each side
- maximum patient weight 330 pounds (150 kg)
- 11 foot power cord with hospital grade plug
- 24 VDC battery charger; 120 VAc, 60 Hz .input



Overall dimensions of the unit:

- width: 26.5 in.
- depth: 38
- height: 44.5

Shipping weight is about 200 pounds.

Supine standing board - \$975



List: \$3,042. MedCon: G. Rifton mobile supine stander.

Features:

- provides standing support and patient transfer aid; allows unrestricted arm

motion

- fold-away hand crank raises patient from horizontal to vertical with any angle in between
- fold-away arm rest
- waist and leg straps
- footboard with lateral support blocks
- four locking casters
- 28 x 45 inch base

Overall dimensions of the unit:

- width: 22 in.
- length: 73 (or standing height)
- height: 32 (reclined)

Shipping weight is about 84 pounds.

Hydrotherapy bath - \$875



List: \$2,400. MedCon: VG. Ille model S hydrotherapy bath. 115 V, 60 Hz.

Features:

- mobile tub with interior dimensions of 28 x 15 x

15.5 inches

- turbine with adjustable height clamp and power on/off toggle switch
- temperature gauge with °C and °F scales
- fill faucet and hose
- 6.5 foot power cord with hospital grade plug

Overall dimensions of the unit:

- width: 16 in.
- length: 32
- height: 38

Shipping weight is about 145 pounds.

Ultrasound therapy unit - \$475



List: \$2,800. MedCon: VG. Enraf Nonius Sonopuls 463 ultrasound therapy unit.

Features:

- pulsed (2.8 ms) or continuous output
- adjustable W/cm²
- adjustable treatment time
- digital display window for min, W/cm², W
- 12 V charger/transformer uses 110 V, 50/60 Hz. current
- to be used with 1 MHz transducer, 5 cm², model no. 1434.80, not included

Overall dimensions of the console:

- width: 13.5 in.
- depth: 10.5
- height: 3.75

Shipping weight is about 15 pounds.

Therapeutic ultrasound - \$575



List: \$1,855. MedCon: VG. Rich-Mar therapeutic ultrasound, 120 V.

Features:

- rolling unit with removable lid, storage cabinet and transport handles
- Rich-Mar C-1 transducer, 1 MHz, 10 cm², beam non-uniformity ratio, 5.5:1
- VI and VI HV B control panels
- power output dial and power gauge; Watts and Watts/cm² scale
- timer
- off/continuous/pulse selector
- polarity selector
- pad connector
- pulse rate selector

- voltage selection dial and DC voltmeter

Overall dimensions of the unit:

- width: 19.5 in.
- depth: 15
- height: 34

Shipping weight is about 60 pounds.

Hydrocollator - \$475



List: \$1,775. MedCon: G. Hydrocollator, no model number, probably made by Chattanooga. 120 V. Good operating condition but missing transport handles and power on indicator light.

Features:

- heating chamber is 20 x 12 x 18 inches
- rack with 4 bars for hanging pads
- power on/off toggle switch
- thermostat adjustment screw
- drain pipe with stopcock
- 6 foot power cord with hospital grade plug


Overall dimensions of the unit:

- width: 24 in.
- depth: 16
- height: 32.5

Shipping weight is about 80 pounds.

Chair scale - \$325 each

List: \$855. MedCon: VG. Detecto weighing chair.



Features:

- 350 lb. capacity
- 16 inch wide seat
- retractable arms and foot rest
- rear wheels and transport handle

Overall dimensions of the unit:

- width: 22 in.
- depth: 27
- height: 41

Shipping weight is about 75 pounds. 2 units in stock.

Adjustable height elevating patient tub - \$1,275



List: \$9,350. MedCon: VG. Arjo elevating hydro massage bath.

Features:

- tub fill valve
- height adjust, up /down
- hot and cold spray hoses
- water temperature gauge
- all stainless steel construction

Shipping weight is about 300 pounds.

Standing aid chair, electrically operated - \$750



List: \$1,900. MedCon: VG. Electrically operated chair, that converts from an armchair to a standing aid.

Features:

- swing away, feeding tray
- power leg rest up/down
- power back rest up/down
- power chair elevation up/down



Shipping weight is about 200 pounds. 2 units in stock.

Sequential compression system - \$375



List: \$1,420. MedCon: VG. Bio Compression Systems, Inc. Sequential Circulator 2004. 120 V, 50/60 Hz.

Features:

- two ports for pressure pads
- pressure increase knob
- pressure gauge, 0-125 mm Hg
- power on/off rocker switch with orange indicator light
- 9 foot power cord with two prong plug

Overall dimensions of the unit:

- width: 12 in.
- depth: 8
- height: 5.25

Shipping weight is about 15 pounds.

Interferential stimulator - \$975



List: \$3,950. MedCon: VG. Chattanooga Intellect 775MP interferential stimulator. 120 V, 60 Hz.

Features:

- designed for treatment of acute and chronic pain, range of motion problems, muscle spasm and circulatory disorders
- beat frequency, 1-200
- red/black electrode balance dial with bar graph LED display
- amplitude adjustment dial
- LED amplitude digital display with red/black output toggle
- 4 sweep modes; constant, 1 set, 80-150 Hz. and 1-10 Hz.
- premodulated, quadpolar, quadpolar with scan, and continuous 5000 Hz modes
- dial-in timer with LED digital display of remaining time upon start
- start pushbutton with indicator light
- power on/off rocker switch
- bumpers/cord wrap on back

Overall dimensions of the unit:

- width: 17.5 in.
- depth: 14.25
- height: 6.5

Shipping weight is about 30 pounds.

Traction machine - \$875



List: \$2,945. MedCon: VG. Dynatronics Citadel traction machine. 120 V, 50/60 Hz.

Features:

- top-mounted traction pulley with cord
- variable speed dial
- dial-in hold and rest traction forces with red LED digital display
- hold, pull rest, and release status indicators (green LED)
- warning indicators (red LED)
- start and reset buttons
- patient remote
- alarms if stopped by remote
- power on/off rocker switch with green indicator light
- bottom plate has drilled holes for mounting unit to traction table

Overall dimensions of the unit:

- width: 10.75 in.
- depth: 14
- height: 10.75

Shipping weight is about 40 pounds.

Heat therapy pump - \$375 each



List: \$1,200. MedCon: VG. Gaymar model TP-500 heat therapy pump. 120 V, 60 Hz., 200 W.

Features:

- inlet and outlet for connecting tubing to therapy pad
- temperature set screw and gauge with °C and °F scales
- fill level indicator
- overtemp warning light
- power on/off rocker switch with green indicator light
- carrying handle
- 8 foot power cord with hospital grade plug; cord tie
- quick instructions and warnings printed on top of pump

Overall dimensions of the unit:

- width: 8 in.
- depth: 5.75

- height: 6.5

Shipping weight is about 20 pounds. 9 units in stock.

Water therapy tank - \$750



List: \$1,900. MedCon: VG. Ille Hydro massage therapy tank

Features:

- approximately 100 gallon capacity
- interior dimensions, 22.5 x 65.5 x 18.5 inches
- heavy gauge stainless steel
- 1/2 hp turbine operates on 115 V, 60 Hz.
- plumbing fixtures included as indicated in the photograph
- drain in bottom of tank and overflow drain

Overall dimensions of the unit:

- width: 23.5 in.
- length: 72
- height: 23 (41 to top of turbine)

Shipping weight is about 230 pounds.

Extremities hydrotherapy tub - \$675



List: \$1,692. MedCon: G. Hydrotherapy physical therapy tub, for extremities, made by Whitehall. For treatment of arms, legs, ankles and knees.

Features:

- built-in pump
- water temperature indicator
- adjustable height support

Interior dimensions:

- length: 25"
- width: 13"
- depth: 15"

Shipping weight is about 140 pounds.

Physical therapy tub - \$975



List: \$1,692. MedCon: G. Hydrotherapy physical therapy tub, for extremities, made by Whitehall.

Features:

- built-in pump
- water temperature indicator
- adjustable height support

Interior dimensions:

- length: 32"
- width: 15"
- depth: 25"

Shipping weight is about 175 pounds.

Traction unit - \$1,375



List: \$4,400. MedCon: VG. Chattanooga TX-1 traction unit. 120 V, 60 Hz.

Features:

- mobile unit with front-locking casters
- digital display of pounds of traction and time remaining

- 0-200 pounds of traction
- treatment modes include preset, progressive intermittent, intermittent, static and progressive static
- progressive modes have up to nine steps
- treatment timer, 1-99 minutes
- selectable traction and rest cycle times
- cord pole has 3 holes for defined heights and a tightening knob to hold at any desired height (to 70 inches)
- cable with patient stop button
- 9 foot power cord with hospital grade plug

Overall dimensions of the unit:

- width: 24.5 in.
- depth: 30
- height: 43 (minimum)

Shipping weight is about 140 pounds.

Hydrocollator - \$475



List: \$1,100. MedCon: G. Ferno ILLE Inferno model 916 moist heat therapy unit. 115 V, 60 Hz.

Features:

- on mobile stand
- variable thermostat controls heating element to 168 degrees F (tested to 55°C)
- temperature adjust screw
- pack rack holds 12 moist heat packs (not included)
- 15 gallon water capacity
- drain valve
- improvised lid was added to replace damaged; gasket doesn't form complete seal
- power on/off rocker switch with indicator light
- 6.5 foot power cord with hospital grade plug

Overall dimensions of the unit:

- width: 27.25 in.

- depth: 16
- height: 31

Shipping weight is about 75 pounds.

Sequential compression pump - \$575



List: \$545. MedCon: VG. Chattanooga Group Pressing Gradient Sequential 3 compression unit (also called Multi 3 or model 4331 or 43000). 110 V, 60 Hz.

Features:

- 72 sec. inflation cycles
- 18 sec. deflation cycles
- pressure control knob, 30-100 mm Hg pressure
- digital display window shows pressure
- two ea., model 43036 three-chambered half leg garment, 20 in., 51 cm; nylon reinforced polyurethane
- two garments may be operated by twin 3-tube ports
- power on/off rocker switch
- carrying handle
- instruction manual included



Overall dimensions of the unit:

- width: 6.75 in.
- depth: 6.25
- height: 5.5

Shipping weight is about 20 pounds.

Patient bathing table - \$775



List: \$2,300. MedCon: VG. Arjo patient bathing table.

Features:

- hydraulic petal raises table height from 21.5 to 34.5 inches
- maximum weight, 330 lbs (150 kg)
- patient transfer pad
- slideaway side rails aid in patient transfer
- drain tilt lever
- drain hole in vinyl
- 4 locking casters

Overall dimensions of the unit:

- width: 72.5 in.
- depth: 27
- height: 30

Shipping weight is about 200 pounds.

Current generator - \$2,275



List: \$7,900. MedCon: VG. Physio Technology OmniStim medium frequency current generator, model 3020. 110/220 V, 60/50 Hz.

Features:

- 4 channel output
- output 200V-P-P Sine curve, 2.5/5 kHz.
- adjustable ramp and on and off times
- digital LED display of time, rate and current; channel selector dial for rate and current
- mode and function dials (see photo)
- dial in timer and start button
- hand control jack
- tabletop or rack mount

Overall dimensions of the unit:

- width: 19 in.
- depth: 13.5
- height: 7.5

Shipping weight is about 25 pounds.

Hydrotherapy unit - \$1,475



List: \$4,850. MedCon: VG. Ferno Ille Hi-Lo Jr. hydrotherapy unit, model 403. 120 V, 60 Hz.

Features:

- two water jets agitate entire bath; 1/2 hp pump
- tub height raises from 20 to 32 inches; convenient for treatment of upper and lower limbs
- LCD digital temperature
- pushbutton start
- adjustable jet force and air content
- disinfect cycle
- drain hose

Overall dimensions of the unit:

- width: 23 in.
- depth: 39.5
- height: 39

Shipping weight is about 200 pounds.

Rehabilitation evaluation system (work simulator) - \$6,000

List: \$38,300. MedCon: VG. BTE work simulator, model 101A. 115 V, 50/60 Hz. This is a hand and upper extremity evaluation system for patient testing with specific tools related to workplace tasks.















The original model was introduced in 1982 and BTE is still in business.

Features:

- automatic graph parameters help assess rehabilitation

outcomes

- maximum strength testing
- consistency of effort
- endurance
- maximum and repetitive lift capacity
- work capacity testing
- FCE components
- control console provides guided sequence of set up operations and prints patient report
- control console has intensity adjustment knob, numerical entry pad and LED displays for distance, time and exercise level
- drum height is adjustable from 37 to 54 inches and tilts for appropriate patient orientations with different tools
- the following tools are included - -
 - BTE 131, small wheel
 - BTE 141, large wheel
 - BTE 162, power grip and finger pinch
 - BTE 122, two-handled reciprocating crank
 - BTE 171, cylindrical grasp
 - BTE 181, multiple handle crossbar
 - BTE 802, large crank handle
 - BTE 191B, extension pole
 - BTE 701 medium crank handle with hand grasp
 - BTE 901, one or two handled grasp (shovel handle)
 - BTE 191 pulley and sand bag
 - BTE 102, one-inch diameter knurled knob
 - BTE 301, two-inch diameter knurled knob
 - BTE 302, three-inch diameter knurled knob
 - BTE 400, small wheel with short crank handle
 - BTE 501, small screwdriver handle
 - BTE 502, medium screwdriver handle
 - BTE 504, large screwdriver handle
 - BTE 601, cylindrical handle
 - BTE 001, range of motion limiter
- 10 foot power cord with hospital grade plug
- Operator's manual included

<p>control console (showing prompts)</p>			
<p>tools included (see list above; proper uses are described in the Operator's Manual)</p>			
<p>Operator's Manual</p>			
			

Overall dimensions of the unit:

- width: 38 in.
- depth: 36
- height: 47

Shipping weight is about 425 pounds.

Bicycle ergometer controller - \$750

List: \$3,200. MedCon: VG. Collins Pedal mate. 115 V, 60 Hz.

Features:



controller for bicycle ergometer compatible with the 5/6 pin cable (see photo)

- manual and heart rate control modes
- heart rate preset, Watt load and timer dials
- analog gauges for actual heart rate, Watt load and pedal rpm

- power on/off dial with indicator lamp
- recorder output and earphone jack

Overall dimensions of the controller:

- width: 16 in.
- depth: 7
- height: 10

Overall dimensions of the controller on stand:

- base: 24 x 24 in.
- height: 40

Shipping weight is about 40 pounds.

Ergometer - \$1,475



List: \$\$\$3,600. MedCon: F. Siemens model 62 03 186 E022E. This appears to be an ergometer that attaches to the end of a gurney or physical therapy table, or attaches over the top of an x-ray or other imaging table.

Features:

- metal pedals with straps attached to enclosed flywheel
- adjustable width; attaches to rails up to 30 in. apart
- 9.5 foot cable with 15 pin connector
- rpm meter is broken and needs to be replaced; connects by standard jack

Overall dimensions of the unit:

- width: 28 in.
- depth: 27

- height: 34

Shipping weight is about 100 pounds.

Rehabilitation exercising system - \$1,750



List: \$58,400. MedCon: F. Loredan Biomedical Lido Active MJ isokinetic rehabilitation system. 120 V, 50/60 Hz. MJ stands for multi-joint. We could not get all of the software to load; may need software update.

Features:

- Lido 486 computer, keyboard, LIDO Impression 3 monitor (13.5 in. diagonal) and Epson FX-870 printer mounted on mobile operator's station with locking casters
- operator's station doubles as a storage board and shelving for test attachments
- dynamometer on rail allows variable positioning of test attachments to patient limbs
- Operations Manual, Service Manual and LIDO Active and LIDO Lift software disks included (see note in List line)



front view shows dynamometer on rail





controller I/O panel
dynamometer cable to J168 is included



various strength and range-of-motion testing attachments



various strength and range-of-motion testing attachments

	Operations and Service Manuals
	LIDO Active and LIDO Lift software

Overall dimensions of the exercising system:

- width: 48 in.
- length: 79
- height: 38.5

Overall dimensions of the operator's station:

- width: 37 in.
- depth: 28
- height: 42

Shipping weight is about 1,100 pounds.

Stretcher scale - \$975 each



List: \$2,650. MedCon: VG. Scale-Tronix model 2001 Sling Scale with sling.

Features:

- double pendulum swivel; twin hoops with hooks will accept a four point stretcher, 72 x 30 in.
- plasticized stretcher included
- shifter handle spreads legs from 17.5 to 36 in. for stability while lifting
- hand-operated hydraulic pump with slow release valve
- transport handle
- rolls on four casters
- clip and rest for securing rolled stretcher
- bed scale with LED digital weight display; reads to 0.1 kg or lb; zero dial and read button

- operating instructions printed on scale casing
- battery charger operates from 115-120 VAC
- 450 lb capacity



Overall dimensions of the unit:

- width: 25.5 in.
- depth: 41
- height: 51.5

Shipping weight is about 140 pounds. 5 units in stock.


Patient lifter with scale - \$1,375



List: \$5,625. MedCon: VG. Detecto Weighmobile model IB-600

Features:

- heavy duty in-bed scale with a capacity of 600 lbs
- twin arms will accept a four point sling; arm raises from 38.5 to 73 in.
- shifter handle spreads legs for stability while lifting
- hand-operated hydraulic pump
- transport handle
- rolls on four casters
- scale with gross/net and lb/kg toggles, zero, tare and measurement lock functions and LCD digital weight display
- 0.2 lb, 0.1 kg sensitivity
- auto power off and auto sleep modes
- Operating Instructions card



Overall dimensions of the unit:

- width: 31.5 in.
- depth: 56
- height: 60

Shipping weight is about 160 pounds.

Stretcher scale - \$975



List: \$4,050. MedCon: VG. Healthometer patient lifter with model 2000 scale.

Features:

- several add-ons to a basic Hoyer lift
- twin hoops with hooks will accept a four point stretcher, 70.5 x 29.5 in.
- shifter handle spreads legs from 19.5 to 32.5 in. stability while lifting
- hand-operated hydraulic pump
- transport handle
- clip and rest for securing rolled stretcher
- scale with lb/kg buttons, tare, reweigh, average lock and LED digital weight display; reads to 0.1 kg or lb
- 12 V battery charger operates from 120 VAC, 60 Hz.
- 400 lb capacity
- rolls on four casters, one locking
- decal with stretcher attachment instructions fixed to main upright tube

Overall dimensions of the unit:

- width: 24.5 in.
- depth: 41
- height: 51.5

Shipping weight is about 140 pounds.

Stretcher scale - \$975



List: \$4,075. MedCon: VG. Healthometer model 2000 KLS patient lift with ProPlus scale.

Features:

- twin arms will accept a four point sling
- shifter handle spreads legs from 18.5 to 35 in. stability while lifting
- hand-operated hydraulic pump
- transport handle
- rolls on four casters
- scale with lb/kg toggle, tare and LED digital weight display; reads to 0.1 kg or lb
- 12 V battery charger operates from 120 VAC, 60 Hz.
- 400 lb capacity

Overall dimensions of the unit:

- width: 23 in.
- depth: 41
- height: 51.5

Shipping weight is about 145 pounds.

Electromechanical patient lifter - \$900



List: \$2,100. MedCon: VG. Linak Hoyer patient lift, model listed as Standard Electrical.

Features:

- Linak J1CB-001 actuator control box and remote control provides slow start and stop for patient safety
- emergency stop button

- LCD window shows battery charge status
- Linak J1BA-001 battery, 24 VDC
- leg spreader controlled by lever with rubber grip
- transport handles
- rolls on four casters
- built-in foot plate and shin support
- 4 hooks for attachment of slings (see photo for orientation)
- shown with large transport sling (included)

Overall dimensions of the unit:

- width: 24 in.
- depth: 48
- height: 45

Shipping weight is about 140 pounds.

Electromechanical patient lifter - \$950 each



List: \$2,600. MedCon: VG. Linak Hoyer patient lift.

Features:

- Linak J2CB-004 actuator control box and remote control provides slow start and stop for patient safety
- Linak J1BA-001 battery, 24 VDC
- green indicator lights for on and charge
- leg spreader controlled by remote control
- transport handles
- rolls on four casters
- 4-point cradle accepts a variety of slings; shown with small sling (included)

Overall dimensions of the unit:

- width: 26 in.
- depth: 46
- height: 54

Shipping weight is about 140 pounds. 2 units in stock.

Hoyer lift - \$475 each



List: \$1,725. MedCon: VG. Hoyer lifter.

Features:

- two point cradle to accept Hoyer slings
- shifter handle spreads legs from 18.5 to 34.5 in. for stability while

lifting

- hand-operated hydraulic pump
- transport handle
- rolls on four casters, two locking

Overall dimensions of the unit:

- width: 24 in.
- depth: 41
- height: 53

Shipping weight is about 70 pounds. 3 units in stock.

Electromechanical patient lifter with scale - \$975



List: \$4,175. MedCon: VG. Medi-Man Rehabilitation Products Medi-Lifter model 622000 patient lift.

Features:

- four point cradle holds sitting position sling
- shifter handle spreads legs from 21 to 40 in. stability while lifting
- electromechanical lift with remote control
- transport handles
- rolls on four casters
- Medi-Man 1000 scale with zero, weigh and lb/kg toggle gives digital weight in LCD window
- battery condition indicator
- uses four TR-10-6 6V rechargeable batteries



Overall dimensions of the unit:

- width: 24 in.
- depth: 43
- height: 52

Shipping weight is about 100 pounds.

Ultrasound therapy unit - \$875



List: \$2,800. MedCon: VG. Enraf Nonius Sonopuls 463 ultrasound therapy unit.

Features:

- 1 MHz transducer, 5 cm², model no. 1434.801
- adjustable W/cm²
- digital display window for min, W/cm², W
- 12 V charger/transformer uses 110 V, 50/60 Hz. current

Overall dimensions of the console:

- width: 13.5 in.
- depth: 10.5
- height: 3.75

Overall dimensions of the case:

- width: 18 in.
- depth: 15.25
- height: 4.25

Shipping weight is about 20 pounds.

Air float system - \$275 each



List: \$670. MedCon: G. Lapidus air float mattress system, made by Hamilton American. Prevents pressure points from building for bed-ridden patients.

Features:

- one/two pads selector switch
- on/off switch
- carrying handle



Shipping weight is about 30 pound. 8 units in stock.

Universal pole clamp - \$30 each



List: \$55 each. MedCon: N. 3 way IV pole/ surgical attachment clamp. Adapts for either IV pole or flat surface mountings. "Quicktrans" clamp, cat # 8C100A, made for IVAC by Pryor products. Black anodized aluminum, machined from 1 solid block.

About 30 units in stock.

Hydrotherapy bath - \$1,475

List: \$5,675. MedCon: VG. Whitehall Lo-Boy 90 gallon hydrotherapy bath. 115 V, 60 Hz.



Features:

- stainless steel tub, 60 x 24 x 18 in. high
- tank rim reinforced with concealed stainless steel rod
- Whitehall model JO-315-SDO turbine with

adjustable positioning device

- dual scale temperature gauge with long stem fitting into retainer
- over-the rim fill spout not included
- 8 foot power cord with hospital grade plug

Overall dimensions of the unit:

- width: 69 in.
- depth: 25
- height: 42

Shipping weight is about 275 pounds.

Hydrotherapy bath, extremity - \$1,275



List: \$4,965. MedCon: VG. Whitehall Hi-Boy 40 gallon hydrotherapy bath. 115 V, 60 Hz.

Features:

- stainless steel tub, 31.5 x 15 x 25 in. high
- tank rim reinforced with concealed stainless steel rod
- Whitehall model JO-75A-SDO turbine with adjustable positioning device
- bottom drain
- overflow drain
- Fahrenheit temperature gauge with long stem
- over-the rim fill spout ready to plumb
- 8.5 foot power cord with hospital grade plug

Overall dimensions of the unit:with him

- width: 36 in.
- depth: 18

- height: 51

Shipping weight is about 150 pounds.

Hydrotherapy bath, extremity - \$975



List: \$5,106. MedCon: VG. Ferno hydrotherapy bath. 115 V, 60 Hz.

Features:

- stainless steel tub, 27.5 x 15 x 18 in. high
- model 100 SAMC turbine, adjustable height, with rubber-covered power on/off switch
- drain plumbing is below tub level (see photo)

Overall dimensions of the unit:

- width: 31 in.
- depth: 16
- height: 46

Shipping weight is about 150 pounds.

Exercise bike - \$215



List: \$549. MedCon: VG. Schwinn DX900 exercise bike.

Features:

- foam wrapped adjustable handlebars
- 31 pound flywheel
- enclosed chain
- LCD displays time, speed, distance, rpm and calories

Overall dimensions of the unit:

- width: 44 in.
- depth: 21.5
- height: 43.5

Shipping weight is about 105 pounds.

Large extremity dry heat fluid therapy - \$1,975



List: \$6,600. MedCon: VG. Henley Fluidotherapy model 210P. 115 V, 60 Hz. Uses fine solid particles in controlled airflow to heat limb for physical therapy purposes. Superficial heat modalities have been reported to be more effective than diathermy and ultrasound in raising muscle temperature. See photograph for dispersion of CELLEX by airflow.

Features:

- leg insertion port (see above photograph) for therapy of knee, ankle or foot
- large capacity for heating and air dispersion of up to 60 lbs of CELLEX media
- adjustable airflow and temperature knobs with heating indicator lamp
- adjustable timer
- 10 foot power cord with hospital grade plug

Overall dimensions of the unit:

- width: 21 in4.
- depth: 44
- height: 35.25

Shipping weight is about 200 pounds.

Upper body ergometer - \$2,750

List: \$12,300. MedCon: VG. Cybex MET-300 arm cranking ergometer. 115 V, 60 Hz.

Features:



- isokinetic and constant power resistance modes
- custom program storage
- adjustable arm crank lengths and seat positions
- 12.5 in. color monitor can display TV or VCR signal with overlaid exercise data; speaker in front of rider
- wheels in front for moving



available programs

Overall dimensions of the unit:

- width: 28 in.
- depth: 61
- height: 60

Shipping weight is about 300 pounds.

Portable battery powered patient lift - \$975



List: \$5,600. MedCon: VG. Ferno Eas-Ille-Lift II portable battery operated patient lift.

Features:

- hand remote with up/down buttons raises seat height from 14.5 to 46 in.
- adjustable side supports/arm rests
- 12 V battery is charged with 120 VAC power cord
- indicator lights for charge needed, power on and battery charging
- rear locking casters

Overall dimensions of the unit:

- width: 32 in.
- depth: 32
- height: 56

Shipping weight is about 200 pounds.

Wheelchair scale - \$750



List: \$2,900. MedCon: VG. Ohaus/Acme Medical Scale Co. model 7401 low profile wheelchair scale.

Features:

- lb/kg toggle to display 0.0 lb or 0.00 kg
- hold button retains measured weight for 2 min.
- uses 105 VAC 50/60 Hz. power to charge six AA nickel-cadmium rechargeable batteries
- roll-on ramps are 5.75 in. wide and 30.75 in. long
- instruction manual included

Overall dimensions of the unit:

- width: 27 in.
- depth: 43
- height: 52

Shipping weight is about 65 pounds.

Full body electric patient lift - \$1,275 each



List: \$4,200. MedCon: G. Andermac Anderlift Plus. 110 V, 50/60 Hz.

This product can be used for patient transport to bed or equipment using a sling sheet and without removing patient from sheet.

Features:

- hand-held hydraulic lift control
- folds to compact storage shape
- maximum capacity, 400 lb
- emergency stop button
- battery charge indicator
- locking casters
- cinch straps but not patient sheet included

Overall dimensions of the unit:

- width: 84 in.
- depth: 24 storage
- height: 43

Shipping weight is about 140 pounds. 2 units in stock.

Sequential compression device - \$375 each



List: \$4,230. MedCon: VG. Kendall model 5325 lymphatic sequential compression device. 115 V, 60 Hz. Cuff and compression hose not included.

Features:

- update of the model 5320 to allow setting of desired pressure by touch pad rather than by dial
- pressure set shown in digital display window
- cycle monitor with ankle, calf and thigh modes
- indicator lamps for power on, sleeve cooling and audible alarms
- recessed carrying handle
- laminated card with operating instructions

Overall dimensions of the unit:

- width: 13 in.
- depth: 8.5
- height: 6.5

Shipping weight is about 20 pounds. 3 units in stock.

Anti-embolism system - \$375 each



List: \$1,925. MedCon: G. Baxter PAS II pulsatile anti-embolism system. 120 V, 60 Hz.

Overall dimensions of the unit:

- width: 12.5 in.
- depth: 6
- height: 6.5

Shipping weight is about 25 pounds. 12 units in stock

Patient lift with scale - \$675

List: \$3,425. MedCon: G. Patient Hoyer lift, with built in. Used to weigh patients who are immobile, or otherwise cannot be moved off of a bed or table. 400 pounds capacity. hydraulically operated lift, through hand pump as shown.



Shipping weight is about 200 pounds.

Physical therapy tub - \$800

List: \$2380. MedCon: VG. All stainless steel physical therapy tub, made by Ille electric manufacturing.



Features:

- adjustable height whirlpool motor
- temperature gauge
- all stainless steel construction

Dimensions:

- length: 52
- width: 21
- depth: 39

Shipping weight is about 200 pounds.

Wheelchair accessible physical therapy tub - \$3,750



List: \$14,210. MedCon: EX. wheelchair accessible physical therapy tub/whirlpool, made by Apollo bath.

Features:

- door at end of tub opens, then seals
- built-in whirlpool motor, with on off switch
- tub drain open/close lever
- built-in water temperature thermometer
- separate hot and cold water inlets and valves
- hand held water sprayer
- in tub seats

- in tub leg rest



Shipping weight is about 400 pounds.

Physical therapy table - \$675



List: \$2,000. MedCon: Demo. Maple top physical therapy table. Top is made of maple wood. Adjusts in height from a minimum of 27 to a maximum of 38 inches in height.

Dimensions:

- length: 66"
- width: 48"
- height: 27-38", variable

Shipping weight is about 300 pounds.

Physical therapy table - \$475

List: \$2,000. MedCon: Demo. Maple top physical therapy table. Top is made of maple wood. Adjusts in height from a minimum of 26 to a maximum of 35 inches in height.



Dimensions:

- length: 60"
- width: 28"
- height: 26-35", variable

Shipping weight is about 200 pounds.

Muscle exerciser - \$375

List: \$1,125. MedCon: G. Muscle exerciser, made by Activeaid. Uses a series of weights and pulleys for exercising muscles. Unit collapses for storage when not in use.



About 6 feet tall overall. Shipping weight 150 pounds.

Lift with patient scale - \$975

List: \$4,350. MedCon: G. Scaletronic model 2001 scale on a Hoyer lift. Capacity is 400 pounds. used to weigh patients either for weight, or while doing a patient transfer. On a Hoyer lift with adjustable legs and hydraulic, hand operated lift.



Shipping weight is about 250 pounds.

Hoyer patient lift - \$475 each



List: \$1,870. MedCon: VG. Hoyer patient lift. For lifting paraplegic patient into and out of whirlpool baths, beds, etc.
Features:

- each from 2 to 6" high
- does not include sling
- adjustable width legs
- 400 pounds capacity
- disassembles for shipping purposes

Shipping weight is about 150 pounds, disassembled. This unit is for anyone who has to administer care for a paraplegic or quadriplegic patient. The long legs of the base will go under either a hospital style beds or a physical therapy whirlpool.

4 units in stock.

Phleboregraph - \$975

List: \$8,200. MedCon: G. Grass Phleboregraph model PRG, 6 channel recorder with accessories, manuals, ink and paper. For diagnosing deep venous thrombosis of lower extremities.



Includes multiple blood pressure cuffs (one for each channel).

Patient trolley - \$700

List: \$2,850. MedCon: VG. Arjo patient trolley. This lift adjusts from a full reclining position to various angles of back and leg support. raises and lowers by a hydraulic foot pump.

Shipped without mattress pad.



Hoyer patient scale / lift - \$875



List: \$4,352. MedCon: VG. Hydraulic patient lift, serial # HL-36028. With digital scale, selectable Kg/ LB, by SRI instruments. 4 point patient bar for suspension of patient sling (no sling). Adjustable spreading legs on base. On casters. Battery operated, with built-in plug-in recharger.

10 lifts of various configurations in stock.

Production Engineering - Medical Equipment Division

6035 E. 38th Avenue, Denver, Colorado 80207 USA (303) 393-7800 FAX (303) 393-1482

© Copyright 1993 - 2011 by Production Engineering. PEMED® is a registered trademark of Production Engineering. MedCon is also a trademark of Production Engineering. All rights reserved. Contact PEMED for license requests. The usual copyright notices and disclaimers. Except for search engines or link listing, no part of this website, including derived information, may be republished without express written permission and payment of fees.

6035 E. 38th Avenue
Denver, Colorado 80207
USA
telephone (303) 393-7800
FAX (303) 393-1482
<http://www.pemed.com>
email: markz@pemed.com

© Copyright 1993 - 2011 by Production Engineering. PEMED® is a registered trademark of Production Engineering. MedCon is also a trademark of Production Engineering. All rights reserved.

6035 E. 38th Avenue
Denver, Colorado 80207
USA
telephone (303) 393-7800
FAX (303) 393-1482
www.pemed.com
markz@pemed.com